+ Evidence in focus

Literature review

Smith-Nephew

Negative pressure wound therapy with instillation and dwell time (NPWTi-d): a review of the evidence

Evidence summary

A systematic literature review of publications about use of NPWTi-d was conducted¹ and the results of high-level (Level 1) studies²⁻⁶ and expert opinion publications⁷⁻¹¹ were summarised. Key findings were:

93% lower-quality studies (116 of 125 are Levels 3 to 5)¹

No statistically significant differences in clinical outcomes (p>0.05) versus NPWT in Level 1 studies²⁻⁵ Expert reviews and guidelines

nighlight a lack of comparative data and inconsistent use of NPWTi-d⁷⁻¹¹

Overview and aim

- Instillation NPWT (NPWTi-d) is a modification of traditional NPWT for the treatment of acute and chronic wound infections¹²
 - It involves instillation of saline, an antiseptic or antibiotic into a sealed wound¹²
- A substantial evidence base of high-level studies supports the use of NPWT in a wide range of wound types, whereas less is known about the effects of NPWTi-d¹²
- A systematic review of the literature about use of NPWTi-d was conducted with particular focus on randomized controlled trials (RCTs) comparing NPWTi-d with NPWT¹
- This report summarizes the results of the highest quality studies (Level 1) that were identified¹ as well as recommendations from consensus groups and health regulators⁷⁻¹¹

Systematic literature review: results for NPWTi-d

- Of the studies published on NPWTi-d up to July 3, 2020, five Level 1 studies (highest quality; Figure 1) were identified in the systematic literature review¹
- Another 120 studies were identified, which were all lower-quality evidence (Levels 2 to 5; Figure 1)¹



Advanced Wound Management, Smith & Nephew Fort Worth, TX 76109

Smith Nephew

Key findings: Level 1 studies

No significant differences in clinical outcomes²⁻⁶

A brief summary of all five Level 1 studies (RCTs) that were identified as part of the systematic literature review is provided below. In general, there were no significant differences in clinical outcomes in studies comparing NPWTi-d with NPWT.²⁻⁵

Lavery LA, et al. 2020²

- This single-center RCT compared healing outcomes in patients with moderate or severe foot infections requiring incision, drainage and intravenous antibiotics
- Patients (n=150) received either NPWT or NPWTi-d with 0.1% polyhexanide-betaine (PHMB) at 30cc/hr (-125mmHg continuous pressure for both) for 16 weeks
- There were no statistically significant differences in wound healing, time to heal, wound dehiscence, re-infection, leg amputation, or hospital re-admission between the two treatments (p>0.05; Figure 2)

'The addition of irrigation to NPWT did not change clinical outcomes in patients with diabetic foot infections.' No statistically significant differences

Figure 2. Clinical outcomes with NPWTi-d and NPWT²

Kim PJ, et al. 2020³

- This multicenter RCT compared NPWTi-d (PHMB; 20min dwell time) with NPWT in acute and chronic wounds requiring operative debridement (-125mmHg continuous pressure)
 - >75% of wounds were on the lower extremities
- A total of 181 patients were treated for 56 days from initial debridement or until the wound was deemed ready for closure or coverage, whichever occurred first
- Results showed no statistically significant difference between the groups in the primary endpoint of number of required inpatient operating room debridements after initial debridement (p>0.05; Figure 3)
- Time to readiness for wound closure/coverage, proportion of wounds closed, and incidence of wound complications were all similar
- Mean decrease in total bacterial count from time of initial surgical debridement to first dressing change was significantly greater with NPWTi-d than NPWT (-0.18 vs 0.60 log₁₀ CFU/g; p=0.02)

1.1 mean number of inpatient debridements

'[Although an unexpected finding] These endpoint results suggest that there was not a strong association between a mean reduction in bioburden and a mean reduction in debridement for the wounds included in this study."

Figure 3. Mean number of inpatient debridements after initial operating room debridement with NPWTi-d and NPWT³

Advanced Wound Management, Smith & Nephew Fort Worth, TX 76109

27167-en V1 1020. Published September 2020. ©2021 Smith+Nephew. [©]Trademark of Smith+Nephew. All Trademarks acknowledged. RNCE21-27041-1220

Smith-Nephew

Davis KE, et al. 2019⁴

- Another single-center RCT compared NPWTi-d using saline (15mL/hr) with two NPWT systems in patients with moderate and severe infected foot wounds (n=30 per group; –125mmHg continuous pressure)
- There were no statistically significant differences between the treatments in any of the assessed clinical outcomes within 12 weeks (p>0.05; Figure 4)
 - Proportion of healed wounds
 - Surgical wound closure rate
 - Number of surgeries
 - Hospital length of stay
 - Time to wound healing

'Our result did not show any benefit when simultaneous irrigation was provided with NPWT as compared to NPWT alone.' Similar results for NPWTi-d and NPWT in all assessed outcomes

Figure 4. Differences in clinical outcomes with NPWTi-d and NPWT⁴

Yang C, et al. 2017⁵

- This RCT compared NPWTi-d (using sodium hypochlorite solution; 10min dwell time) with NPWT (-125 mmHg continuous pressure) in 20 chronic wounds (10 per group) after one week of treatment:
 - Compared with baseline, there was a reduction in biofilm-protected bacteria in the NPWTi-d group (p<0.05) and a slight increase in the NPWT group (p=0.46)
 - However, when the results for each treatment group were compared, the difference was not significant (p=0.11; Figure 5)
 - Use of NPWTi-d did not significantly reduce planktonic bacteria concentrations after initial debridement (p=0.16)
 - Change in wound size was similar with both treatments (Figure 5)

NPWTi-d significantly reduced biofilm-protected bacteria versus baseline (p<0.05), but not planktonic bacteria; similar changes in wound size with both treatments

> "...it remains unclear as to what degree of biofilm-protected bacteria reduction must be achieved to elicit a clinical result."

Figure 5. Changes in wound size and reductions in biofilm-protected bacteria with NPWTi-d and NPWT $^{\rm 5}$

Kim PJ, et al. 20156

- Rather than comparing with NPWT, this RCT evaluated the effects of different NPWTi-d solutions on treatment outcomes in patients (n=100) with wounds that required hospital admission and surgical debridement
 - Most wounds were located on the foot or lower leg (>70%)
- Use of saline was shown to be as effective as antiseptic (PHMB; 20min dwell time) but had shorter time to final surgical procedure than antiseptic (p=0.038)
- There were no statistically significant differences between saline and antiseptic for number of operating room visits, length of hospital stay, and number of wounds that closed and remained closed at 30-day follow-up (p≥0.05)

'Perhaps the positive effect observed with negative-pressure wound therapy with instillation reflects this intermitent negative-pressure application phenomenon rather than anything to do with instillation of a solution.'

Saline as effective

Figure 6. Outcomes using saline or antiseptic as the instillation solution with NPWTi-d^6 $\,$

Advanced Wound Management, Smith & Nephew Fort Worth, TX 76109

Smith_{Nephew}

What experts say

- · Several clinical guidelines and expert reviews advocate that large Level 1 RCTs are needed to evaluate the efficacy of NPWTi-d versus NPWT in the clinical setting⁷⁻¹¹
- The most recent consensus guidelines by Kim PJ and colleagues acknowledge that more comparative data are needed to determine the effectiveness of NPWTi-d⁸
- One health regulator in the UK has produced specific advice on NPWTi-d⁹
 - It states that although NPWTi-d appears more effective than moist wound care and NPWT for acutely infected or chronic non-healing wounds, it is difficult to draw clear conclusions due to inconsistent use and a lack of RCTs⁹

Summary

- Results of Level 1 studies comparing NPWT show no significant differences in clinical outcomes between the two treatments²⁻⁵
- Most studies evaluating NPWTi-d are lower-quality evidence (93% of 125 studies are Levels 3 to 5)¹
- high-level (Level 1) evidence supporting the use of NPWTi-d¹

For detailed product information, including indications for use, contraindications, precautions and warnings, please consult the product's applicable Instructions for Use (IFU) prior to use.

References

1. Ebohon S. A systematic literature review of instillation therapy. Smith+Nephew Evidence analysis report. EO/AWM/RENASYS/002/v3. 27 July 2020. 2. Lavery LA, Davis KE, La Fontaine J, et al. Does negative pressure wound therapy with irrigation improve clinical outcomes? A randomized clinical trial in patients with diabetic foot infections. Am J Surg. 27 Feb 2020 [Epub ahead of print]. 3. Kim PJ, Lavery LA, Galiano RD, et al. The impact of negative-pressure wound therapy with instillation on wounds requiring operative debridement: Pilot randomised, controlled trial. Int Wound J. 21 Jun 2020. [Epub ahead of print]. 4. Davis KE, La Fontaine J, Farrar D, et al. Randomized clinical study to compare negative pressure wound therapy with simultaneous saline irrigation and traditional negative pressure wound therapy for complex foot infections. Wound Repair Regen. 2020;28(1):97–104. 5. Yang C, Goss SG, Alcantara S, Schultz G, Lantis li JC. Effect of negative pressure wound therapy with instillation on bioburden in chronically infected wounds. 2017;29(8):240-246. 6. Kim PJ, Attinger CE, Oliver N, et al. Comparison of outcomes for normal saline and an antiseptic solution for negative-pressure wound therapy with instillation. Plast Reconstr Surg. 2015;136(5):657e-664e. 7. McKanna M, Geraci J, Hall K, et al. Clinician panel recommendations for use of negative pressure wound therapy with instillation. Ostomy Wound Manage. 2016;62(4):S1-S14. 8. Kim PJ, Attinger CE, Constantine T, et al. Negative pressure wound therapy with instillation: International consensus guidelines update. Int Wound J. 2020;17(1):174–186. 9. National Institute for Health and Care Excellence. Medtech Innovation Briefing 189. The V.A.C. Veraflo Therapy system for infected wounds. 2 Sep 2019. Available at: https://www.nice.org.uk/advice/mib189. Accessed 05 Aug 2020. 10. Gupta S, Gabriel A, Lantis J, Téot L. Clinical recommendations and practical guide for negative pressure wound therapy with instillation. Int Wound J. 2016;13(2):159–174. 11. Anghel EL, Kim PJ. Negative-pressure wound therapy: A comprehensive review of the evidence. Plast Reconstr Surg. 2016; 138(3 Suppl):1295–375. (Abstract only). 12. Apelqvist J, Willy C, Fagerdahl AM, et al. Negative Pressure Wound Therapy – overview, challenges and perspectives. J Wound Care. 2017;26(Suppl 3):S1-S113.